


1.9011.6.22 / 1.9011.T.33

Foam Rubber Insert for Blood Tests



Foam rubber insert 1.9011.6.22 (left) and 1.9011.T.33 (left) – Insert 1.9011.T.33 with carrier 1.9350.T.1B and test tubes (right)

Purpose: Protecting the inserted test tubes

Material: Pu foam rubber (washing temperature: 120°C)

Different dimensions available:

Foam rubber cap for inserting in carriers 1.9224 / DN110:

Has capacity for: 16 tubes each insert
One carrier can accommodate two inserts
with 32 tubes in total

Dimensions: 83 mm × 115 mm
Weight: 0.0132 kg
Art.no.: 1.9011.6.22

Foam rubber cap for inserting in carriers 1.9350.T.xxH:

Has capacity for: 9 tubes each side of
an insert
One insert can accommodate 18 tubes
in total

Dimensions: 123 mm × 100 mm
Weight: 0.0246 kg
Art.no.: 1.9011.T.33

We reserve the right of changes!